

Cologuard® bee na'alkaahí bił Ach'íidíil biyi'di łóód binaałniih beinít'ijh bá ná'náł'ijh

Cologuard bee na'alkaahí Medicare dóó aláahdi béeso ach'áqah nideinilí lq'í da **yik'éda'asti'**.

\$0 Diné nídabidi'nél'ijhí béeso bich'áqah nidaa'níilí lq'í da Cologuard bee na'alkaahí bá **\$0** yik'é nidayiiléh.^{1*}

96% Kéyah nitsaa bił hahoodzo t'áa attsohjił, Cologuard bee diné nídabidi'nél'ijhí **96%** bilááh ná'náł'ijh ats'aazisdóó bąqah ál'íní nijiléhígíi yik'é nidayiiléh.^{1*}



Ła' bínááhwiidoo'áatgo kojji' baa nanítáh **Cologuard.com/insurance** doodago kojji' **1-844-870-8870 hólne'**.

Béeso nich'áqah naa'nít bił haz'ání nida'ídikiid t'áa ákódaat'éhígíi ná nabídiilkiid bee níka'iilyeedgo bíghah.

Béeso ách'áqah nideinil bił nidahaz'ání lq'í da Cologuard ná'nél'ijh t'áa aníi át'éeé lq'í hodoo'niid **bikéédóó ach'íidíil bii' al'thgo naalkaah bee níná'nél'ijhí yik'éda'asti'**.^{2,3†}




Haash Yił'éego Chojooł'ijh Naalkidígíi: Cologuard.com/use

Cologuard bee na'alkaahí éi haghandi chqo'ijigo doo nanitł'ah da


Bik'ehgo Ada'doolníígíi Attso Hazhó'ó Baa Hane'ígíi éi Sits'aa' Cologuard biyi' si'ánígíi biyi'.

1 DAH NÍDAAH




- K'ad bíghah ninz'ijgo Cologuard choidíit'ijit.
- Łigai bee na'at'áhi bikáá' dah ná'nídahí bikáá' dah díi'áat.
- Bikáá' dah ná'nídaahí bikáá' dah nídaah.

2 BÍNÍGHAAS




- Łigai bee na'at'áhi biyi'jił adíichijit.
- Tsin áts'óózigíi bee nichaan (chqo') bidiighas dóó bee si'ánidéé' attsoozígíi biyi' góne' abídíiyit.

3 NÍDÍILŁOH Í BIDÁ'DIDÍI'ÁÁL



- Tó t'igai bee na'at'áhi bee ha'didíit'ijit.
- Bee na'at'áhi nitsaaígíi bii' shínichqo'ígíi biih díi'áat.

4 BÍ'DÍILŁAH Í BIŁ Í'DÍILÍIL



- Bídeiljahí t'áa átał attsoł hadíidíilíit.
- Díi bee si'ání Exact Sciences UPS® bee ínídíidíilíit.



Bee í'inilyaaígíi t'áa jiik'eh bił índíidíilíit.


1. UPS deididoo'áatgo bee bił nihwiní'aah kojji' **Cologuard.com/UPS** doodago kojji' **1-844-870-8870 hólne'**.
2. Bee si'ání UPS bich'ij' díi'áat. T'áa ákóqoh bił haz'ánígíi dóó da'ılkeed góne' nehajáhígíi nił bééhodooz'ijitgo kohjił **Cologuard.com/UPS** baa nanítáh.

*Atk'idáqá' diné nábidi'nél'ijhí bee bąqah háájilígíi bik'ehgo Exact Sciences aniid shq'shin ha'niih.

†Diné nídabidi'nél'ijhí lq'í da bá, Yas Niłt'ees 1, 2023 yootkáat góne' bee adiníish, dóó t'áa sahdii bee hada'dít'éhí naaltsoos bee ókeedí bee lq'í azl'ígíi bitígíi doo bideet'i' da.^{2,3}

Ni Cologuard® bee na'alkaahí bee béédahoozinígíí

Exact Sciences bee í'inilaaígíí néidii'á bikéédóó nika'análwo'í
ats'íís yaa'áhalyání bee ná béédahoozinígíí neidiiniít 8-10 yiskáqjít'
hólqo dooleet.



Bee na'alkaahí Cologuard éi nichaan DNA biyi' **lahgo ádzaaígíí dóó diit**
(tóód binaatniih bééhozinígíí) héidées'íí'⁴

- Naanii ni'ídzoh bee bééhoozingo t'áa ha'nínígíí Cologuard bee na'alkaahí tóód binaatniih
bee éédahoziní doo bik'í'neetqá da. Ach'íidíil biyi'di tóód binaatniih doo nee hólqo da.⁴
- + Atnáh ni'ídzoh bee bééhoozingo t'áa ha'nínígíí Cologuard bee na'alkaahí tóód binaatniih bee
éédahoziní bik'í'neetqá. Atnáh ni'ídzoh bee bééhozinígíí tóód binaatniih nee hólqo **doo** t'áa
ha'níi da.⁴Yah anáa'ádaáh biniiyé nika'análwo'í ats'íís yaa'áhalyání bit nihwiní'aah.

Diné ta' doo ákót'éeégóo (doo ákót'ée da) Atnáh ni'ídzoh dóó naanii ni'ídzoh bee bá béédahoozih.
Náa'ákóne' náánidi'nél'ínígíí baqah nika'análwo'í ats'íís yaa'áhalyání átahjít' nabídítkid.⁴

Bee baa dahane'í ta' náána nit bééhodooziit dóó
nidaalkidí kojít': [Cologuard.com/UnderstandingResults](https://www.cologuard.com/UnderstandingResults)

Binhajj' Éédahozinígíí dóó Bee Anáhodoot'ih Bee Baa Dahane'í Bódahólníhígíí

Cologuard® dóó Cologuard Plus® bee nida'alkaahí éi dahóyáanii eę dóó bilááh béédaahaaígíí hái ach'íidíil
tóód binaatniih bee bich'j' anídahodoot'ih bit nidawólnígíí DNA bee bééhoziní haashjít' daat' éhígíí dóó chqá'
biyi' diit béédahoozihgo bee dabidi'nél'j' biniiyé. Atk'idqá' dah ach'íidíil tóód binaatniih doodago ats'íistoo'
bii' tóód doo t'áa ákót'éeégóo nooséígíí, neezgaigo ajichj' tóód binaatniih bee nik'ih hoolzhiizh nit'ée'go,
atk'idqá' nighan haz'áqgi ach'íidíil tóód binaatniih doodago binahjít' bee éédahoziní bit oolchítígíí, doodago
náána ta' ach'íidíil tóód binaatniih bee ná'nél'j' at'éhígíí éi na'alkaah bee ná'nél'j' yaa ch'óhwiní'ánígíí
bee nanidiská'ígíí t'áa aníí at'éego bee bééhoozingo Cologuard ádaalyaaígíí t'áadoo choinít'íní. Díí bee
ádahoot'éhígíí háidíil shjít' nidadeet'i'go nika'análwo'í ats'íís yaa'áhalyání bee bich'j' hadíidzih.

Cologuard bee na'alkaahí bee éédahoozinígíí t'áa baa hasti'go bee it hane' dooleet. Bee na'alkaahí atnáh
ni'ídzoh bee bééhozinígíí tóód binaatniih hólónígíí doo yee hoot'aah da. Diné nidabidi'nél'j'j'íí bee na'alkaah
atnáh ni'ídzoh bee bá béédahoozinígíí ach'íidíil bii' alt'ihgo bee na'alkaah bich'j' bida'iilniih dooleet. Bee
na'alkaahí naanii ni'ídzoh bee bééhozinígíí tóód binaatniih ádinígíí doo yee hoot'aah da. Diné nidabidi'nél'j'j'íí
bee na'alkaah naanii ni'ídzoh bee bá béédahoozinígíí hahgo nináadabi'dilkaahígíí biniiyé be' azéé' íit'íní yee
yít ahít' halne' dooleet. Atnáh ni'ídzoh dóó naanii ni'ídzoh doo ákót'éeégóo béédahoozihígíí anéego bíghah.
Rx t'éiyá.

Binhajj' Éédahoziní: 1. Internal Data on File. Exact Sciences. 2. FAQs about Affordable Care Act Implementation Part 51, Families First Coronavirus Response Act and Coronavirus Aid, Relief, and Economic Security Act Implementation. US Department of Health & Human Services. Published January 10, 2022. Accessed March 30, 2026. <https://www.hhs.gov/guidance/document/faqs-about-affordable-care-act-implementation-part-51>. 3. CMS. Calendar year (CY) 2023 Medicare physician fee schedule final rule. Published November 1, 2022. Accessed March 30, 2026. <https://www.cms.gov/newsroom/fact-sheets/calendar-year-cy-2023-medicare-physician-fee-schedule-final-rule>. 4. Exact Sciences. 2026. Cologuard Patient Brochure.

EXACT SCIENCES CORPORATION | 5505 Endeavor Lane, Madison, WI 53719
ExactSciences.com | ExactLabs.com | 1-844-870-8870

Cologuard dóó Exact Sciences éi Exact Sciences Corporation bibee ééhóziní bit
ha'dilyaaígíí át'é. Náána ta' bibee éédahoziní t'áa atsoh daabi daniinígíí daa bí.
© 2026 Exact Sciences Corporation. Bee adahazt'i'í t'áa atsoh baa hóchxj'.

M-US-CG-07088_LL_April 2026

 **cologuard**[®]
by EXACT SCIENCES