

Kuaj ntsuas txog kab mob khees xaws rau plab hnyuv nrog rau qhov kev sim kuaj Cologuard®

Kev sim kuaj mob Cologuard **tau txais kev pab them nqi duav roos** los ntawm Medicare thiab cov kev tuav pov hwm uas tseem ceeb tshaj plaws.

\$0

Cov neeg mob uas muaj is saws las tuav pov hwm feem coob lawv them \$0 rau kev sim kuaj mob Cologuard^{1*}

96%

Thooob plaws hauv teb chaws, ntau tshaj **96%** ntawm cov neeg mob Cologuard tsis tau them tus nqi uas koj rho nyiaj hauv hnab tshos them (out-of-pocket cost) rau qhov kev kuaj ntsuas mob^{1*}



Kawm paub ntau ntxiv nyob rau ntawm Cologuard.com/insurance los sis hu rau **1-844-870-8870**. Peb muaj peev xwm pab koj nug koj lub tuam txhab is saws laws tuav pov hwm uas siv cov lus nug kom raug.

Cov kws is saws las tuav pov hwm feem coob lawv **pab them nqi duav roos rau ib qho kev siv teeb tsum tshuaj xyuas hnyuv laus taug qab xyuas ntxiv mus** tom qab ua tiav ib qho kev kuaj ntsuas mob uas siv kev kuaj ntsuas mob Cologuard uas pom tias muaj mob huam tuaj.^{2,3†}




Kev sim kuaj mob Cologuard nws **yooj yim** siv nyob rau ntawm vaj tsev

Cov Lus Qhia Meej Tseeb Txog Txhua Yam yog muaj nyob rau hauv Lub Npov Cologuard

qhov vis dis aus qhia txog
Kev Yuav Siv Tau Li Cas:
Cologuard.com/use

1


ZAUM



- Siv qhov kev sim kuaj mob Cologuard thaum twg koj npaj txhij lawm.
- Muab koj lub taub ntim qauv piv txwv tsos dawb tso rau ntawm lub qhov viv tso quav.
- Zaum rau ntawm lub qhov viv tso quav.

2


DAUS



- Tso quav rau hauv lub taub ntim qauv piv txwv tsos dawb ntawd.
- Muab daim ntawv nplaum daus koj qhov quav (cov quav) thiab muab tso rov qab rau hauv lub me nyuam taub ntim qauv uas muab los ntawm qhov cuab yeej pab cawm kev raug mob ntawd.

3


TSAU DEJ THIAB FOOB



- Muab tshuaj pov thaiv kom tsis txhob lwj tsw hom ua kua tso rau hauv lub taub ntim qauv piv txwv ntawd.
- Muaj koj qhov qauv piv txwv kaw cia rau hauv lub taub ntim qauv uas loj.

4

LO DAIM NTAWV NPLAUM MUAB LUS QHIA RAU THIAB XA MUS



- Sau lus qhia rau ntawm cov ntawv nplaum muab lus qhia ob qho tib si.
- Muab qhov cuab yeej pab cawm kev raug mob no xa rov qab mus rau Exact Sciences los ntawm UPS®.



Xa koj qhov qauv piv txwv rov qab yam tsis them nqi dab tsi li


- Teem ib lub sij hawm mus nqa ntawm UPS nyob rau ntawm Cologuard.com/UPS los sis hu rau **1-844-870-8870**
- Nqa koj qhov cuab yeej pab cawm kev raug mob mus rau UPS. Tshawb nrhiav cov khw muag khoom hauv lub zej zos thiab cov xuab moos uas koj mus tso rau tau nyob rau ntawm Cologuard.com/UPS

*Exact Sciences qhov kev kwv yees tsis ntev los nod yog saib raws li kev sau nqi los ntawm tus neeg mob uas sau tseg rau hauv ntau ntawv.

†Rau cov neeg mob feem coob, qhov no yuav pib siv tau nyob rau Lub Ib Hlis Ntuj 1, 2023, thiab siv tsis tau rau cov phiaj xwm tuav pov hwm fab kev lag luam uas muaj qhov txheej xwm kev zam.^{2,3}

Koj qhov kev kuaj mob Cologuard® uas tau los ntawm kev sim kuaj mob

Koj tus kws muab kev saib xyuas kho mob rau fab kev noj qab haus huv yuav tsum muaj koj cov kev kuaj tau los nyob rau 2 lub lim tiam tom qab Exact Sciences txais tau koj qhov qauv piv txwv lawm.



Qhov kev sim kuaj mob Cologuard tshawb nrhiav txog **DNA** thiab **ntshav** uas **pauv hloov** (cov tsos mob ntawm tus kab mob khees xaws) hauv koj cov quav⁴



Ib qhov kev kuaj tau los uas pom tias tsis muaj mob txhais hais tias qhov kev sim kuaj mob Cologuard nrhiav tsis pom muaj cov tsos mob ntawm tus kab mob khees xaws uas muaj feem tshwm sim tau. Ntxim li koj yuav tsis muaj tus kab mob khees xaws (cancer) rau txoj hnyuv laus.⁴



Ib qhov kev kuaj tau los uas pom tias tus mob huam tuaj txhais hais tias qhov kev sim kuaj mob Cologuard nrhiav pom muaj cov tsos mob ntawm tus kab mob khees xaws uas muaj feem tshwm sim tau. Ib qhov kev kuaj tau los uas pom tias tus mob huam tuaj ntawd **tsis tau** txhais hais tias koj muaj tus kab mob khees xaws.⁴ Teem sij hawm nrog koj tus kws muab kev saib xyuas kho mob rau fab kev noj qab haus huv rau ib qho kev soj qab taug lw xyuas mus ntxiv.

Qee cov neeg tau txais kev kuaj tau los ntawm kev sim kuaj mob uas pom tias tus mob huam loj tuaj thiab tsis muaj mob yuam kev (tsis raug). Nco ntsoov nug koj tus kws muab kev saib xyuas kho mob rau fab kev noj qab haus huv txog cov kauj ruam ntxiv mus rau kev kuaj ntsuas mob.⁴

Tshawb nrhiav lus qhia paub thiab cov vis dis aus ntau ntxiv tau nyob rau ntawm: **Cologuard.com/UnderstandingResults**

Lus Qhia thiab Cov Ntaub Ntawv Qhia Txog Kev Phom Sij Tseem Ceeb

Qhov kev sim kuaj mob Cologuard yog npaj los mus kuaj ntsuas mob rau cov neeg laus hnuv nyoog 45 xyoo thiab laus dua ntawd uas nyob rau kev phom sij hauv sim txog rau kev mob tus kab mob khees xaws (cancer) rau txoj hnyuv los ntawm kev kuaj txhom qee cov cim ntsig txog DNA thiab cov ntshav uas nyob rau hauv cov quav. Tsis txhob siv yog hais tias koj muaj adenomas, muaj kab mob hnyuv uas ua rau hnyuv liab thiab o thiab qee cov kev mob ib tiam dhau rau tiam, los sis ib qho keeb kwm hais txog kab mob khees xaws (cancer) rau txoj hnyuv ntawm tus kheej los sis tsev neej. Ib qho kev sim kuaj mob Cologuard tsis yog ib qho kev hloov pauv rau kev siv teeb tsum tshuaj xyuas hnyuv laus hauv cov neeg mob uas muaj kev phom sij siab. Kev sim kuaj mob Cologuard hauv cov neeg laus hnuv nyoog 45-49 xyoo yuav raug kwv yees saib raws li qhov kev kawm tshawb fawb fab kev kho mob loj txog ntawm cov neeg mob hnuv nyoog 50 xyoo thiab laus dua ntawd.

Qhov kev kuaj tau los ntawm kev sim kuaj mob Cologuard yuav tsum raug muab tshab ntxhais lub ntsiab yam xyuam xim. Ib qhov kev kuaj tau los ntawm kev kuaj ntsuas mob pom tias muaj tus mob huam tuaj lawm ntawd yeej tseem qhia tsis tau laib qhov tseeb tias yog mob khees xaws lawm. Yuav tsum tau muab cov neeg mob uas muaj ib qhov kev kuaj tau los ntawm kev kuaj ntsuas mob pom tias muaj tus mob huam tuaj lawm ntawd xa tauj mus ua ib qho kev siv teeb tsum tshuaj xyuas hnyuv laus. Ib qhov kev kuaj tau los ntawm kev kuaj ntsuas mob pom tias tsis muaj mob ntawd yeej tseem qhia tsis tau laib qhov tseeb tias tsis muaj mob khees xaws. Cov neeg mob uas muaj ib qhov kev kuaj tau los ntawm kev kuaj ntsuas mob pom tias tsis muaj mob yuav tsum tau sib tham nrog lawv tus kws kho mob thaum lawv xav tau kev sim kuaj ntsuas mob ib zaug ntxiv. Cov kev kuaj tau los ntawm kev sim kuaj ntsuas mob pom tias muaj tus mob huam tuaj lawm uas tsis tseeb thiab cov kev kuaj tau los ntawm kev sim kuaj ntsuas mob pom tias tsis muaj mob uas tsis tseeb muaj peev xwm tshwm sim tau. Nyob rau hauv ib qho kev kawm tshawb fawb txog fab kev kho mob, muaj txog 13% ntawm cov neeg uas tsis muaj tus kab mob khees xaws tau txais ib qhov kev kuaj tau los ntawm kev sim kuaj ntsuas mob pom tias muaj tus mob huam tuaj lawm (tsis tseeb) thiab 8% ntawm cov neeg uas muaj tus kab mob khees xaws tau txais ib qhov kev kuaj tau los ntawm kev sim kuaj ntsuas mob pom tias tsis muaj mob (tsis tseeb). Rx nkaus xwb

Cov Ntaub Ntawv Ua Pov Thawj: 1. Internal Data on File. Exact Sciences. 2. FAQs about Affordable Care Act Implementation Part 51, Families First Coronavirus Response Act and Coronavirus Aid, Relief, and Economic Security Act Implementation. US Department of Health & Human Services. Luam Tawm Lub Ib Hlis Ntuj 10, 2022. Nkag Cuag Tau Lub Tsib Hlis Ntuj 2, 2024. <https://www.hhs.gov/guidance/document/faqs-about-affordable-care-act-implementation-part-51>. 3. CMS. Calendar year (CY) 2023 Medicare physician fee schedule final rule. Luam Tawm Lub Kaum Ib Hlis Ntuj 1, 2022. Nkag Cuag Tau Lub Tsib Hlis Ntuj 2, 2024. <https://www.cms.gov/newsroom/fact-sheets/calendar-year-cy-2023-medicare-physicianfee-schedule-final-rule>. 4. Exact Sciences. 2022. Phau Ntawv Qhia Txog Tus Neeg Mob Ntawm Cologuard.

EXACT SCIENCES CORPORATION | 5505 Endeavor Lane, Madison, WI 53719
ExactSciences.com | ExactLabs.com | 1-844-870-8870

Cologuard thiab Exact Sciences yog ib cov cim kev lag luam uas tau txais kev sau npe raug cai ntawm Exact Sciences Corporation.

Tag nrho lwm cov cim lag luam puav leej yog peev txheej ntawm lawv cov tswv ntawd li ntiag tug.

© 2024 Exact Sciences Corporation. Txwv tsis pub yuam txhua txoj cai. MS-US-CG-05477-LL-Hmong-White-August 2024

 **cologuard®**
by EXACT SCIENCES

**EXACT
SCIENCES**